





Introduction to Real-Time PCR and CFX System



Charlene Lu
2021/02/2`
Bio-Rad Laboratories, Taiwan




Generations of PCR

1st Generation

2nd Generation


3rd Generation



Gel Electrophoresis
(qualitative)

Real-Time PCR
(indirect quantification)


Droplet Digital PCR
(absolute quantification)



Polymerase Chain Reaction (PCR)

Targeted DNA replication using thermostable DNA polymerase

DENATURE






Polymerase Chain Reaction (PCR)


The use of two primers allows targeting of specific sequences

Primers are complementary to opposite strands of target region but not complementary to any other sequences

ANNEAL PRIMERS




50 – 65°C



Polymerase Chain Reaction (PCR)


Taq polymerase synthesizes DNA complementary to template in 5' to 3' direction

EXTEND STRANDS 72°C





Polymerase Chain Reaction (PCR)



Each cycle of PCR doubles the number of progeny DNA duplexes (which can then act as template as well)


1 cycle = 2^1 copies of starting template

25 cycles = 2^{25} copies of starting template
(~33554432)

Amplification is exponential, but the exponential increase is limited:


- 1 A linear increase follows exponential
 - Eventually plateaus

Real-Time PCR allows us to 'see' the exponential phase so we can calculate how much we started with.




PCR – 1st Generation



- Amplify target DNA with end point analysis to distinguish products
 - No relationship between end point and starting target copies



End point analysis




qPCR Principle





- Record target-specific fluorescence during the PCR and then determination C_q Value for sample
- An inverse, linear relationship exists between the logarithm of the initial target quantity and the C_q
- Signal threshold determines C_q value (formerly C_t)
- Logarithmic base 2 scale

Determining Threshold and Cq

$$X_n = X_0(1+E)^n$$



• Normal view versus Log view



C_T value v.s. concentration

- 1 cycle = 2 fold difference
- 3.32 cycles \cong 10 fold difference
- Assumes 100% efficiency

$$Y = N_0 2^n$$




$$Y = N_0 (1+E)^n$$



Quantitative PCR (qPCR)

Create a **standard curve** with 10-fold serial dilutions of PCR product – assign arbitrary values

Compare values from standards with values for unknown sample



Standard Curves




A standard curve allows you to determine the efficiency of your assay

C_T is linear with the log of starting copy number (standard curve)

Correlation Coefficient: 0.99 Slope: -3.396 Intercept: 28.373 $Y = -3.396 X + 28.373$
PCR Efficiency: 97.0%

□ Unknowns
○ Standards



The **slope** of the standard curve can be directly correlated to the **efficiency** of the reactions:

r = is a measure of how well the actual data fit to the standard curve.
= (explained variation/total variation)


Aim for R value (Correlation Coefficient) of ≥ 0.98

Efficiency (E) = $[10^{(-1/\text{slope})}] - 1$
when slope = -3.32, Efficiency = 100%

Efficiency Requirements

- Assay design
 - Secondary structures?
 - Primer-dimers?
 - Nonspecific amplification
- Pipetting?
- Dilution errors?

- Assay design
 - Secondary structures?
 - Primer-dimers?
 - Nonspecific amplification
- Pipetting?
- Dilution errors?
- PCR inhibitors?



E	Slope
0.5	-5.679
0.6	-4.899
0.7	-4.339
0.8	-3.917
0.9	-3.587
1	-3.322
1.1	-3.103
1.2	-2.920
1.3	-2.765

$$\text{Slope} = - [1/\log (1+E)]$$

$$\log (1+E) = - (1/\text{slope})$$


$$E = 10^{[-1/\text{slope}]} - 1$$

Aim for Efficiency Values:

Good = 90 – 110%



Fantastic = 95 – 105%

- **Threshold Cycle values (C_T) have a direct relationship to the amount of starting template**
- **Check space between C_T values follows correct relationship (100% efficiency) $2^n = \text{fold dilution}$**
- **Efficiency of reactions between 90-110%**
- **R value should be ≥ 0.98**




What are the most common detection strategies used for Real-Time PCR?

- The fluorescent molecules can be used
 - Non-specific DNA binding dyes
 - SYBR[®] Green I
 - Specific Hybridization Probes/Primers
 - TaqMan[™]



Real-Time Chemistry: SYBR Green



- SYBR Green I binds dsDNA
- fluorescence increases when bound to dsDNA
- As dsDNA accumulates, more dye binds the DNA and the fluorescence increase



Intercalation dyes: SYBR and EvaGreen

Advantages

- Intercalation dyes-based assays are more economical than probes-based ones
- Allow melt curve analysis to test assay specificity


Disadvantages

- You can not multiplex with dyes
- Dyes do not discriminate and will bind any dsDNA in your sample
 - This includes primer dimers
 - Problems if assay is not specific

Melt Curve Analysis


Fluorescence decrease as the temperature increase:

1. DNA strands start to separate
2. SYBR green loses its binding to the DNA
3. Fluorescence rapidly decreases




Melt Curve Analysis


- The melting temperature of the amplicon can easily be detected.
- Contaminating DNA, primer dimer or false priming is seen as an additional peak.



Melt Curve Analysis




Multiple wells, same assay, single peaks at high melting temperatures ($+80^{\circ}\text{C}$). This means my assay is specific!



Multiple wells, same assay, multiple peaks! This means off-target amplification and/or primer dimers!

Real-Time Chemistry: Probes (TaqMan)

- Target specific hybridization probe
- 5' reporter and 3' quencher
- Utilizes FRET (Fluorescence Resonance Energy Transfer) quenching
- Relies on a validated primer set with high specificity



* heat for BHQs


Dark Quenchers :Ground state or static quenching

Energy is dissipated as heat

Molecular interactions inhibit fluorescence

Mechanism of probe chemistry

1. During PCR, probe hybridizes to target sequence
2. Probe is partially displaced during extension
3. Probe cleaved by 5'-3' nuclease activity of polymerase, releasing reporter from quencher
4. Illuminated reporter exhibits unquenched fluorescence





Hybridization probes


Advantages

- Allows multiplexing by using different reporter dye
- Removes concerns about primer dimers and off-target amplifications (**More specific**)
- SNP genotyping application

Disadvantages

- More expensive than dye-based assays
- To test assay specificity you'll need to run products on gels

Gel confirmation



SYBR vs. Probes: How do I choose?

34

- Is your sample limited or rare?
- Is time-to-results a very significant constraint?
- Are you differentiating similar sequences?
- Do you have specificity concerns?
- Do you need multiple data points out of the same sample?



Sso Advanced Universal SYBR®
Green Supermix
1725270



Sso Advanced Universal Probes
Supermix
1725280

Real-time PCR Sample Preparation

SYBR Green Chemistry

Component	Volume per reaction	Final concentration
iQ SYBR Green Supermix	25 μ l	1X
Primer 1	x μ l	100 nM–500 nM
Primer 2	x μ l	100 nM–500 nM
Sterile water	x μ l	
DNA template	x μ l	
Total Volume	50 μl	


Probe Chemistry

Component	Volume per reaction	Final concentration
iQ Supermix	25 μ l	1X
Primer 1	x μ l	100 nM–500 nM
Primer 2	x μ l	100 nM–500 nM
Probe	x μ l	100 nM–500 nM
Sterile water	x μ l	
DNA template	x μ l	
Total Volume	50 μl	



Hot Start

PCR

Melting Curve



3 Step PCR



3 min Polymerase Activation
 40 min Denaturing
 40 min Annealing
 40 min Extension
 7 min Final Extension


130 min total in incubations

Initial Ramp to 95°
 40 X Ramp down to 60°
 40 X Ramp up to 72°
 39 X Ramp up to 95°

46 min total in ramping

Total time 2 hours 56 minutes

2 Step PCR



3 min Polymerase Activation
7 min Denaturing
40 min Annealing

50 min total in incubations

Initial Ramp to 95°
40 X Ramp down to 60°
39 X Ramp up to 95°

46 min total in ramping

Total time 1 hour 36 minutes



General QPCR Working Process -- Data analysis

- ***Relative Quantity (ΔCT)***

- Not normalized
- Normalization accomplished via equal loading of samples
- Post analysis normalization

- ***Normalized Expression ($\Delta\Delta CT$)***


- Accounts for loading differences
- Usually normalize to reference gene
- Relative quantity of Target is normalized by the relative quantity of the reference genes



Relative Quantity (ΔCT)

	Target
Tissue #1:	22
Tissue #2:	24
<hr/>	
Delta Ct:	$24 - 22 = 2$

$$\text{Fold induction} = 2^2 = 4$$




Normalized Expression ($\Delta\Delta CT$)

	Reference	Target
Tissue #1 (Control):	21	22
Tissue #2 (Test) :	20	24
<hr/>		
1 st Delta	Delta Ct #1:	$22 - 21 = 1$
	Delta Ct #2:	$24 - 20 = 4$
<hr/>		
2 nd Delta	Delta Ct:	$1 - 4 = -3$


Fold induction = $2^{-(-3)} = 8$

Relative Quantification

- ❖ Problem with the $\Delta\Delta CT$
- ❖ Slopes are not parallel




Pfaffl modification




	Primer set #1 Reference	Primer set #2 Target
Tissue #1 (Control):	21	22
Tissue #2 (Test) :	20	24

(From Standard curve)	Efficiency:	90% = 1.9	100% = 2
	Delta Ct:	20-21 = -1	24-22 = 2

$$\text{Fold induction} = \frac{2_{\text{target}}^{\text{deltaCt}_{\text{target}} (24-22 = 2)}}{1.9_{\text{reference}}^{\text{deltaCt}_{\text{reference}} (20-21 = -1)}} = \frac{4}{0.56} = 7.1$$



- Δ Ct method: (*no reference gene*)
 - Fold induction : 4
- $\Delta\Delta$ Ct method: (*reference gene*)
 - Fold induction : 8
- Pfaffl modification: (*reference gene and efficiency*)
 - Fold induction : 7.1



Remember...



Real-Time PCR is not 'cookbook chemistry' - a real-time instrument will not optimize your experiments for you

However, once you do **optimize your reactions**, you will get **reproducible, accurate results**

- Primer design considerations
 - Assay type (SYBR vs Probe/TaqMan)
 - Exonic vs intron spanning
 - Selection of primary genome viewer (NCBI, Ensembl, UCSC)
 - SNP locations
 - Selection of best fragments for optimal transcript coverage and maximum intron size
 - Secondary structures (UNAFold)
 - Specificity (BiSearch)


Choosing between intron/exon coverage

Exonic designs

- Assay completely located within a single exon
- Amplifies both cDNA & gDNA
- gDNA can be used as positive control
- Can be applied to all genes

Intron-spanning designs

- At least one intron between forward & reverse primer
- Limits co-amplification of contaminating gDNA
- Not possible for single exon genes or genes with small introns (~ 9% of the genome)
- May not help for genes with processed pseudogenes (~15% of genome)
- Decreased design space



Primer3Plus Primer Design

- <http://www.primer3plus.com>

Primer3Plus
pick primers from a DNA sequence

[More...](#) [Source Code](#)
[Help](#) [About](#)

Load server settings: *Select primer pairs to detect the given template sequence. Optionally targets and included/excluded regions can be specified.*

Task:


Main | **General Settings** | **Advanced Settings** | **Internal Oligo** | **Penalty Weights** | **Advanced Sequence**

Sequence Id:

[Paste template sequence below](#) Or upload sequence file:

- PCR settings
 - mix_salt_divalent
 - mix_salt_monovalent
 - mix_dntp_conc
 - mix_dna_conc
 - mix_pcr_temperature
- amplicon properties
 - amplicon_size
- Tm prediction
 - SantaLucia 1998
- primer properties
 - primer_size
 - primer_tm
 - primer_max_tm_diff
 - primer_gc
 - primer_max_3prime_gc
 - primer_max_3_poly
 - primer_max_self_end
 - primer_max_self_any
 - primer_pair_max_compl_end

Specificity Analysis



BLAST® Basic Local Alignment Search Tool

Home Recent Results Saved Strategies Help

▸ NCBI BLAST Home

BLAST finds regions of similarity between biological sequences. [more...](#)

New DELTA-BLAST, a more sensitive protein-protein search

BLAST Assembled RefSeq Genomes


Choose a species genome to search, or [list all genomic BLAST databases](#).

- [Human](#)
- [Mouse](#)
- [Rat](#)
- [Arabidopsis thaliana](#)
- [Oryza sativa](#)
- [Bos taurus](#)
- [Danio rerio](#)
- [Drosophila melanogaster](#)
- [Gallus gallus](#)
- [Pan troglodytes](#)
- [Microbes](#)
- [Apis mellifera](#)

Basic BLAST

Choose a BLAST program to run.

nucleotide blast	Search a nucleotide database using a nucleotide query <i>Algorithms:</i> blastn, megablast, discontinuous megablast
protein blast	Search protein database using a protein query <i>Algorithms:</i> blastp, psi-blast, phi-blast, delta-blast
blastx	Search protein database using a translated nucleotide query
tblastn	Search translated nucleotide database using a protein query
tblastx	Search translated nucleotide database using a translated nucleotide query



Assays – What are your choices?

- Design new assays
 - Take advantage of the latest information and design tools and do the proper wet-lab validation in order to meet MIQE guidelines
- Use existing custom designed assays
 - Possible unknown design criteria
 - Have they been validated to MIQE standards?
- Use pre-designed assays from a vendor
 - Unknown design criteria (for most)
 - In-silico validation vs. wet-lab validation
 - **Not all assays are created equally!**
 - Bio-Rad's PrimePCR assays/plate arrays
 - **www.Bio-Rad.com/PrimePCR**

PrimePCR Assays

Available Species

Human	Rabbit
Mouse	Rhesus Monkey
Rat	Yeast
Chicken	Zebrafish
Arabidopsis	Cow
Pig	Dog
	Rice

The screenshot shows the Bio-Rad website for PrimePCR assays. The page layout includes a top navigation bar with links for 'Products', 'Applications & Technologies', 'Promotions', 'News & Events', and 'Support'. The main content area is titled 'PrimePCR™ Assays and Panels' and features a 'Why PrimePCR?' section with a 'Get Help' button and a 'My PrimePCR Products' section. A search bar is visible at the bottom of the page.

Also: Human lncRNA assays!

www.bio-rad.com/PrimePCR

PrimePCR Pathway & Collection Panels


Diseases

Cancer (133)
 Cardiac hypertrophy (7)
 Cardiovascular diseases (15)
 Congenital, hereditary, and neonatal diseases and abnormalities (73)
 Cystic fibrosis (37)
 Digestive system diseases (13)
 Endocrine system diseases (20)
 Hemic and lymphatic diseases (51)
 Immune system diseases (53)
 Mental disorders (11)
 Musculoskeletal diseases (6)
 Nervous system diseases (16)
 Nutritional and metabolic diseases (46)
 Respiratory tract diseases (31)
 Skin and connective tissue diseases (19)

96-well	384-well
ABI 7300 - 96	ABI 7900 - 384
ABI 7500 - 96 Fast	ABI QuantStudio - 384
ABI 7500 - 96 Std	ABI ViiA 7 - 384
ABI 7900 - 96 Fast	Bio-Rad CFX384
ABI 7900 - 96 Std	Bio-Rad CFX384 Touch
ABI QuantStudio - 96 Fast	Roche LC480 - 384
ABI QuantStudio - 96 Std	
ABI StepOnePlus - 96	
ABI ViiA 7 - 96 Fast	
ABI ViiA 7 - 96 Std	
Bio-Rad CFX Connect - 96	
Bio-Rad CFX96	
Bio-Rad CFX96 Touch	
Bio-Rad Chromo4 - 96	
Bio-Rad iQ - 96	
Bio-Rad iQ5 - 96	
Bio-Rad MyiQ - 96	
Bio-Rad MyiQ2 - 96	
Bio-Rad Opticon - 96	
Bio-Rad Opticon 2 - 96	
Eppendorf Realplex - 96	
Roche LC480 - 96	
Stratagene Mx - 96	

Metabolism


amino acid metabolism (4)
 carbohydrate metabolism (4)
 lipid metabolism (5)
 nucleotide metabolism (3)
 regulation of lipid metabolism (3)
 steroid metabolism (5)
 vitamin and cofactor metabolism (4)
 xenobiotic metabolism (6)




Additional tips: Experimental Design


- Technical vs. Biological replication
 - Reviewers typically like technical triplicates
 - Smaller fold changes require a higher n of samples
- Multiple reference genes can stabilize basal gene expression
 - Critical if you measure small fold changes

Biological Replicates



RT-qPCR Samples Technical Replicates





Additional tips: Experimental Design

- Volumes matter!
 - St Dev goes up as volume of reaction goes down
 - Avoid pipetting small volumes (1-2 ul)
 - Better to add more of a diluted template
- Errors increase with each pipetting step
- Make up a master mix of everything minus template
 - Usually total samples +10% - always need overage
 - Minimize error, lower total pipetting steps

Example well:

- 10 ul of 2x Supermix
(dNTPs/Taq/Mg⁺⁺/Buffer)
- 1 ul of 20x diluted F+R primer mix
- 3 ul of water
- 5 ul template



Additional tips: Essential Controls

- Non-template control (NTC)
 - Add water in lieu of template
 - Any signal before cycle 38 should be checked
 - Either due to primer dimers or contamination
 - Best to have one NTC well per assay (primer set)
- No reverse transcription (NRT)
 - Checks for DNA carry over in your RNA purification
 - Add equivalent amount of non-converted RNA as you did cDNA
 - ΔCq to cDNA should be above ~ 4 cycles
 - More important on a per sample/prep basis



Additional tips: Essential Controls

- Negative control (Neg)
 - Add template from another species/system/knock-out that does not contain your sequence of interest
 - Differs from NTC – checks for *specificity*
 - Essential to catch off target effects early
- Positive control (Pos)
 - Add your template of interest
 - Look for a signal $<Cq=30$
 - Used primarily in troubleshooting inhibitory samples

Experiment design

Sample maximization method

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk-1 Actb	Unk-1 Actb	Unk-1 Actb	Unk-9 Actb	Unk-9 Actb	Unk-9 Actb	Unk-17 Actb	Unk-17 Actb	Unk-17 Actb	Unk-25 Actb	Unk-25 Actb	Unk-25 Actb
B	Unk-2 Actb	Unk-2 Actb	Unk-2 Actb	Unk-10 Actb	Unk-10 Actb	Unk-10 Actb	Unk-18 Actb	Unk-18 Actb	Unk-18 Actb	Unk-26 Actb	Unk-26 Actb	Unk-26 Actb
C	Unk-3 Actb	Unk-3 Actb	Unk-3 Actb	Unk-11 Actb	Unk-11 Actb	Unk-11 Actb	Unk-19 Actb	Unk-19 Actb	Unk-19 Actb	Unk-27 Actb	Unk-27 Actb	Unk-27 Actb
D	Unk-4 Actb	Unk-4 Actb	Unk-4 Actb	Unk-12 Actb	Unk-12 Actb	Unk-12 Actb	Unk-20 Actb	Unk-20 Actb	Unk-20 Actb	Unk-28 Actb	Unk-28 Actb	Unk-28 Actb
E	Unk-5 Actb	Unk-5 Actb	Unk-5 Actb	Unk-13 Actb	Unk-13 Actb	Unk-13 Actb	Unk-21 Actb	Unk-21 Actb	Unk-21 Actb	Unk-29 Actb	Unk-29 Actb	Unk-29 Actb
F	Unk-6 Actb	Unk-6 Actb	Unk-6 Actb	Unk-14 Actb	Unk-14 Actb	Unk-14 Actb	Unk-22 Actb	Unk-22 Actb	Unk-22 Actb	Unk-30 Actb	Unk-30 Actb	Unk-30 Actb
G	Unk-7 Actb	Unk-7 Actb	Unk-7 Actb	Unk-15 Actb	Unk-15 Actb	Unk-15 Actb	Unk-23 Actb	Unk-23 Actb	Unk-23 Actb	NTC	NTC	NTC
H	Unk-8 Actb	Unk-8 Actb	Unk-8 Actb	Unk-16 Actb	Unk-16 Actb	Unk-16 Actb	Unk-24 Actb	Unk-24 Actb	Unk-24 Actb			

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk-1 GAPDH	Unk-1 GAPDH	Unk-1 GAPDH	Unk-9 GAPDH	Unk-9 GAPDH	Unk-9 GAPDH	Unk-17 GAPDH	Unk-17 GAPDH	Unk-17 GAPDH	Unk-25 GAPDH	Unk-25 GAPDH	Unk-25 GAPDH
B	Unk-2 GAPDH	Unk-2 GAPDH	Unk-2 GAPDH	Unk-10 GAPDH	Unk-10 GAPDH	Unk-10 GAPDH	Unk-18 GAPDH	Unk-18 GAPDH	Unk-18 GAPDH	Unk-26 GAPDH	Unk-26 GAPDH	Unk-26 GAPDH
C	Unk-3 GAPDH	Unk-3 GAPDH	Unk-3 GAPDH	Unk-11 GAPDH	Unk-11 GAPDH	Unk-11 GAPDH	Unk-19 GAPDH	Unk-19 GAPDH	Unk-19 GAPDH	Unk-27 GAPDH	Unk-27 GAPDH	Unk-27 GAPDH
D	Unk-4 GAPDH	Unk-4 GAPDH	Unk-4 GAPDH	Unk-12 GAPDH	Unk-12 GAPDH	Unk-12 GAPDH	Unk-20 GAPDH	Unk-20 GAPDH	Unk-20 GAPDH	Unk-28 GAPDH	Unk-28 GAPDH	Unk-28 GAPDH
E	Unk-5 GAPDH	Unk-5 GAPDH	Unk-5 GAPDH	Unk-13 GAPDH	Unk-13 GAPDH	Unk-13 GAPDH	Unk-21 GAPDH	Unk-21 GAPDH	Unk-21 GAPDH	Unk-29 GAPDH	Unk-29 GAPDH	Unk-29 GAPDH
F	Unk-6 GAPDH	Unk-6 GAPDH	Unk-6 GAPDH	Unk-14 GAPDH	Unk-14 GAPDH	Unk-14 GAPDH	Unk-22 GAPDH	Unk-22 GAPDH	Unk-22 GAPDH	Unk-30 GAPDH	Unk-30 GAPDH	Unk-30 GAPDH
G	Unk-7 GAPDH	Unk-7 GAPDH	Unk-7 GAPDH	Unk-15 GAPDH	Unk-15 GAPDH	Unk-15 GAPDH	Unk-23 GAPDH	Unk-23 GAPDH	Unk-23 GAPDH	NTC	NTC	NTC
H	Unk-8 GAPDH	Unk-8 GAPDH	Unk-8 GAPDH	Unk-16 GAPDH	Unk-16 GAPDH	Unk-16 GAPDH	Unk-24 GAPDH	Unk-24 GAPDH	Unk-24 GAPDH			

Actb

Good DATA

GAPDH

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk-1 P53	Unk-1 P53	Unk-1 P53	Unk-9 P53	Unk-9 P53	Unk-9 P53	Unk-17 P53	Unk-17 P53	Unk-17 P53	Unk-25 P53	Unk-25 P53	Unk-25 P53
B	Unk-2 P53	Unk-2 P53	Unk-2 P53	Unk-10 P53	Unk-10 P53	Unk-10 P53	Unk-18 P53	Unk-18 P53	Unk-18 P53	Unk-26 P53	Unk-26 P53	Unk-26 P53
C	Unk-3 P53	Unk-3 P53	Unk-3 P53	Unk-11 P53	Unk-11 P53	Unk-11 P53	Unk-19 P53	Unk-19 P53	Unk-19 P53	Unk-27 P53	Unk-27 P53	Unk-27 P53
D	Unk-4 P53	Unk-4 P53	Unk-4 P53	Unk-12 P53	Unk-12 P53	Unk-12 P53	Unk-20 P53	Unk-20 P53	Unk-20 P53	Unk-28 P53	Unk-28 P53	Unk-28 P53
E	Unk-5 P53	Unk-5 P53	Unk-5 P53	Unk-13 P53	Unk-13 P53	Unk-13 P53	Unk-21 P53	Unk-21 P53	Unk-21 P53	Unk-29 P53	Unk-29 P53	Unk-29 P53
F	Unk-6 P53	Unk-6 P53	Unk-6 P53	Unk-14 P53	Unk-14 P53	Unk-14 P53	Unk-22 P53	Unk-22 P53	Unk-22 P53	Unk-30 P53	Unk-30 P53	Unk-30 P53
G	Unk-7 P53	Unk-7 P53	Unk-7 P53	Unk-15 P53	Unk-15 P53	Unk-15 P53	Unk-23 P53	Unk-23 P53	Unk-23 P53	NTC	NTC	NTC
H	Unk-8 P53	Unk-8 P53	Unk-8 P53	Unk-16 P53	Unk-16 P53	Unk-16 P53	Unk-24 P53	Unk-24 P53	Unk-24 P53			

P53

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk-1 EGFR	Unk-1 EGFR	Unk-1 EGFR	Unk-9 EGFR	Unk-9 EGFR	Unk-9 EGFR	Unk-17 EGFR	Unk-17 EGFR	Unk-17 EGFR	Unk-25 EGFR	Unk-25 EGFR	Unk-25 EGFR
B	Unk-2 EGFR	Unk-2 EGFR	Unk-2 EGFR	Unk-10 EGFR	Unk-10 EGFR	Unk-10 EGFR	Unk-18 EGFR	Unk-18 EGFR	Unk-18 EGFR	Unk-26 EGFR	Unk-26 EGFR	Unk-26 EGFR
C	Unk-3 EGFR	Unk-3 EGFR	Unk-3 EGFR	Unk-11 EGFR	Unk-11 EGFR	Unk-11 EGFR	Unk-19 EGFR	Unk-19 EGFR	Unk-19 EGFR	Unk-27 EGFR	Unk-27 EGFR	Unk-27 EGFR
D	Unk-4 EGFR	Unk-4 EGFR	Unk-4 EGFR	Unk-12 EGFR	Unk-12 EGFR	Unk-12 EGFR	Unk-20 EGFR	Unk-20 EGFR	Unk-20 EGFR	Unk-28 EGFR	Unk-28 EGFR	Unk-28 EGFR
E	Unk-5 EGFR	Unk-5 EGFR	Unk-5 EGFR	Unk-13 EGFR	Unk-13 EGFR	Unk-13 EGFR	Unk-21 EGFR	Unk-21 EGFR	Unk-21 EGFR	Unk-29 EGFR	Unk-29 EGFR	Unk-29 EGFR
F	Unk-6 EGFR	Unk-6 EGFR	Unk-6 EGFR	Unk-14 EGFR	Unk-14 EGFR	Unk-14 EGFR	Unk-22 EGFR	Unk-22 EGFR	Unk-22 EGFR	Unk-30 EGFR	Unk-30 EGFR	Unk-30 EGFR
G	Unk-7 EGFR	Unk-7 EGFR	Unk-7 EGFR	Unk-15 EGFR	Unk-15 EGFR	Unk-15 EGFR	Unk-23 EGFR	Unk-23 EGFR	Unk-23 EGFR	NTC	NTC	NTC
H	Unk-8 EGFR	Unk-8 EGFR	Unk-8 EGFR	Unk-16 EGFR	Unk-16 EGFR	Unk-16 EGFR	Unk-24 EGFR	Unk-24 EGFR	Unk-24 EGFR			

EGFR

Experiment design

Someone did lay-out like this...

	1	2	3	4	5	6	7	8	9	10	11	12		1	2	3	4	5	6	7	8	9	10	11	12
A	Unk Actb 1	Unk Actb 1	Unk Actb 1	Unk Actb 3	Unk Actb 3	Unk Actb 3	Unk Actb 5	Unk Actb 5	Unk Actb 5	Unk Actb 7	Unk Actb 7	Unk Actb 7	A	Unk Actb 8	Unk Actb 8	Unk Actb 8	Unk Actb 10	Unk Actb 10	Unk Actb 10	Unk Actb 12	Unk Actb 12	Unk Actb 12	Unk Actb 14	Unk Actb 14	
B	Unk GAPDH 1	Unk GAPDH 1	Unk GAPDH 1	Unk GAPDH 3	Unk GAPDH 3	Unk GAPDH 3	Unk GAPDH 5	Unk GAPDH 5	Unk GAPDH 5	Unk GAPDH 7	Unk GAPDH 7	Unk GAPDH 7	B	Unk GAPDH 8	Unk GAPDH 8	Unk GAPDH 8	Unk GAPDH 10	Unk GAPDH 10	Unk GAPDH 10	Unk GAPDH 12	Unk GAPDH 12	Unk GAPDH 12	Unk GAPDH 14	Unk GAPDH 14	Unk GAPDH 14
C	Unk P53 1	Unk P53 1	Unk P53 1	Unk P53 3	Unk P53 3	Unk P53 3	Unk P53 5	Unk P53 5	Unk P53 5	Unk P53 7	Unk P53 7	Unk P53 7	C	Unk P53 8	Unk P53 8	Unk P53 8	Unk P53 10	Unk P53 10	Unk P53 10	Unk P53 12	Unk P53 12	Unk P53 12	Unk P53 14	Unk P53 14	Unk P53 14
D	Unk EGFR 1	Unk EGFR 1	Unk EGFR 1	Unk EGFR 3	Unk EGFR 3	Unk EGFR 3	Unk EGFR 5	Unk EGFR 5	Unk EGFR 5	Unk EGFR 7	Unk EGFR 7	Unk EGFR 7	D	Unk EGFR 8	Unk EGFR 8	Unk EGFR 8	Unk EGFR 10	Unk EGFR 10	Unk EGFR 10	Unk EGFR 12	Unk EGFR 12	Unk EGFR 12	Unk EGFR 14	Unk EGFR 14	Unk EGFR 14
E	Unk Actb 2	Unk Actb 2	Unk Actb 2	Unk Actb 4	Unk Actb 4	Unk Actb 4	Unk Actb 6	Unk Actb 6	Unk Actb 6	Unk Actb 8	Unk Actb 8	Unk Actb 8	E	Unk Actb 11	Unk Actb 11	Unk Actb 11	Unk Actb 13	Unk Actb 13	Unk Actb 13	Unk Actb 15	Unk Actb 15	Unk Actb 15	NTC	NTC	NTC
F	Unk GAPDH 2	Unk GAPDH 2	Unk GAPDH 2	Unk GAPDH 4	Unk GAPDH 4	Unk GAPDH 4	Unk GAPDH 6	Unk GAPDH 6	Unk GAPDH 6	Unk GAPDH 8	Unk GAPDH 8	Unk GAPDH 8	F	Unk GAPDH 11	Unk GAPDH 11	Unk GAPDH 11	Unk GAPDH 13	Unk GAPDH 13	Unk GAPDH 13	Unk GAPDH 15	Unk GAPDH 15	Unk GAPDH 15	NTC	NTC	NTC
G	Unk P53 2	Unk P53 2	Unk P53 2	Unk P53 4	Unk P53 4	Unk P53 4	Unk P53 6	Unk P53 6	Unk P53 6	Unk P53 8	Unk P53 8	Unk P53 8	G	Unk P53 11	Unk P53 11	Unk P53 11	Unk P53 13	Unk P53 13	Unk P53 13	Unk P53 15	Unk P53 15	Unk P53 15	NTC	NTC	NTC
H	Unk EGFR 2	Unk EGFR 2	Unk EGFR 2	Unk EGFR 4	Unk EGFR 4	Unk EGFR 4	Unk EGFR 6	Unk EGFR 6	Unk EGFR 6	Unk EGFR 8	Unk EGFR 8	Unk EGFR 8	H	Unk EGFR 11	Unk EGFR 11	Unk EGFR 11	Unk EGFR 13	Unk EGFR 13	Unk EGFR 13	Unk EGFR 15	Unk EGFR 15	Unk EGFR 15	NTC	NTC	NTC
A	Unk Actb 15	Unk Actb 15	Unk Actb 15	Unk Actb 17	Unk Actb 17	Unk Actb 17	Unk Actb 19	Unk Actb 19	Unk Actb 19	Unk Actb 21	Unk Actb 21	Unk Actb 21	A	Unk Actb 24	Unk Actb 24	Unk Actb 24	Unk Actb 26	Unk Actb 26	Unk Actb 26	Unk Actb 28	Unk Actb 28	Unk Actb 28	Unk Actb 30	Unk Actb 30	
B	Unk GAPDH 15	Unk GAPDH 15	Unk GAPDH 15	Unk GAPDH 17	Unk GAPDH 17	Unk GAPDH 17	Unk GAPDH 19	Unk GAPDH 19	Unk GAPDH 19	Unk GAPDH 21	Unk GAPDH 21	Unk GAPDH 21	B	Unk GAPDH 24	Unk GAPDH 24	Unk GAPDH 24	Unk GAPDH 26	Unk GAPDH 26	Unk GAPDH 26	Unk GAPDH 28	Unk GAPDH 28	Unk GAPDH 28	Unk GAPDH 30	Unk GAPDH 30	
C	Unk P53 15	Unk P53 15	Unk P53 15	Unk P53 17	Unk P53 17	Unk P53 17	Unk P53 19	Unk P53 19	Unk P53 19	Unk P53 21	Unk P53 21	Unk P53 21	C	Unk P53 24	Unk P53 24	Unk P53 24	Unk P53 26	Unk P53 26	Unk P53 26	Unk P53 28	Unk P53 28	Unk P53 28	Unk P53 30	Unk P53 30	
D	Unk EGFR 15	Unk EGFR 15	Unk EGFR 15	Unk EGFR 17	Unk EGFR 17	Unk EGFR 17	Unk EGFR 19	Unk EGFR 19	Unk EGFR 19	Unk EGFR 21	Unk EGFR 21	Unk EGFR 21	D	Unk EGFR 24	Unk EGFR 24	Unk EGFR 24	Unk EGFR 26	Unk EGFR 26	Unk EGFR 26	Unk EGFR 28	Unk EGFR 28	Unk EGFR 28	Unk EGFR 30	Unk EGFR 30	
E	Unk Actb 16	Unk Actb 16	Unk Actb 16	Unk Actb 18	Unk Actb 18	Unk Actb 18	Unk Actb 20	Unk Actb 20	Unk Actb 20	NTC	NTC	NTC	E	Unk Actb 23	Unk Actb 23	Unk Actb 23	Unk Actb 25	Unk Actb 25	Unk Actb 25	Unk Actb 27	Unk Actb 27	Unk Actb 27	NTC	NTC	NTC
F	Unk GAPDH 16	Unk GAPDH 16	Unk GAPDH 16	Unk GAPDH 18	Unk GAPDH 18	Unk GAPDH 18	Unk GAPDH 20	Unk GAPDH 20	Unk GAPDH 20	NTC	NTC	NTC	F	Unk GAPDH 23	Unk GAPDH 23	Unk GAPDH 23	Unk GAPDH 25	Unk GAPDH 25	Unk GAPDH 25	Unk GAPDH 27	Unk GAPDH 27	Unk GAPDH 27	NTC	NTC	NTC
G	Unk P53 16	Unk P53 16	Unk P53 16	Unk P53 18	Unk P53 18	Unk P53 18	Unk P53 20	Unk P53 20	Unk P53 20	NTC	NTC	NTC	G	Unk P53 23	Unk P53 23	Unk P53 23	Unk P53 25	Unk P53 25	Unk P53 25	Unk P53 27	Unk P53 27	Unk P53 27	NTC	NTC	NTC
H	Unk EGFR 16	Unk EGFR 16	Unk EGFR 16	Unk EGFR 18	Unk EGFR 18	Unk EGFR 18	Unk EGFR 20	Unk EGFR 20	Unk EGFR 20	NTC	NTC	NTC	H	Unk EGFR 23	Unk EGFR 23	Unk EGFR 23	Unk EGFR 25	Unk EGFR 25	Unk EGFR 25	Unk EGFR 27	Unk EGFR 27	Unk EGFR 27	NTC	NTC	NTC

BAD DATA

Without inter-run calibration!

Inter-run Calibrator

Gene maximization method

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk Actb 1	Unk Actb 1	Unk Actb 1	Unk Actb 3	Unk Actb 3	Unk Actb 3	Unk Actb 5	Unk Actb 5	Unk Actb 5	Unk Actb 7	Unk Actb 7	Unk Actb 7
B	Unk GAPDH 1	Unk GAPDH 1	Unk GAPDH 1	Unk GAPDH 3	Unk GAPDH 3	Unk GAPDH 3	Unk GAPDH 5	Unk GAPDH 5	Unk GAPDH 5	Unk GAPDH 7	Unk GAPDH 7	Unk GAPDH 7
C	Unk P53 1	Unk P53 1	Unk P53 1	Unk P53 3	Unk P53 3	Unk P53 3	Unk P53 5	Unk P53 5	Unk P53 5	Unk P53 7	Unk P53 7	Unk P53 7
D	Unk EGFR 1	Unk EGFR 1	Unk EGFR 1	Unk EGFR 3	Unk EGFR 3	Unk EGFR 3	Unk EGFR 5	Unk EGFR 5	Unk EGFR 5	Unk EGFR 7	Unk EGFR 7	Unk EGFR 7
E	Unk Actb 2	Unk Actb 2	Unk Actb 2	Unk Actb 4	Unk Actb 4	Unk Actb 4	Unk Actb 6	Unk Actb 6	Unk Actb 6	NTC	NTC	NTC
F	Unk GAPDH 2	Unk GAPDH 2	Unk GAPDH 2	Unk GAPDH 4	Unk GAPDH 4	Unk GAPDH 4	Unk GAPDH 6	Unk GAPDH 6	Unk GAPDH 6	NTC	NTC	NTC
G	Unk P53 2	Unk P53 2	Unk P53 2	Unk P53 4	Unk P53 4	Unk P53 4	Unk P53 6	Unk P53 6	Unk P53 6	NTC	NTC	NTC
H	Unk EGFR 2	Unk EGFR 2	Unk EGFR 2	Unk EGFR 4	Unk EGFR 4	Unk EGFR 4	Unk EGFR 6	Unk EGFR 6	Unk EGFR 6	NTC	NTC	NTC

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk Actb 8	Unk Actb 8	Unk Actb 8	Unk Actb 10	Unk Actb 10	Unk Actb 10	Unk Actb 12	Unk Actb 12	Unk Actb 12	Unk Actb 12	Unk Actb 12	Unk Actb 12
B	Unk GAPDH 8	Unk GAPDH 8	Unk GAPDH 8	Unk GAPDH 10	Unk GAPDH 10	Unk GAPDH 10	Unk GAPDH 12	Unk GAPDH 12	Unk GAPDH 12	Unk GAPDH 12	Unk GAPDH 12	Unk GAPDH 12
C	Unk P53 8	Unk P53 8	Unk P53 8	Unk P53 10	Unk P53 10	Unk P53 10	Unk P53 12	Unk P53 12	Unk P53 12	Unk P53 12	Unk P53 12	Unk P53 12
D	Unk EGFR 8	Unk EGFR 8	Unk EGFR 8	Unk EGFR 10	Unk EGFR 10	Unk EGFR 10	Unk EGFR 12	Unk EGFR 12	Unk EGFR 12	Unk EGFR 12	Unk EGFR 12	Unk EGFR 12
E	Unk Actb 26	Unk Actb 26	Unk Actb 26	Unk Actb 28	Unk Actb 28	Unk Actb 28	Unk Actb 30	Unk Actb 30	Unk Actb 30	NTC	NTC	NTC
B	Unk GAPDH 26	Unk GAPDH 26	Unk GAPDH 26	Unk GAPDH 28	Unk GAPDH 28	Unk GAPDH 28	Unk GAPDH 30	Unk GAPDH 30	Unk GAPDH 30	NTC	NTC	NTC
C	Unk P53 26	Unk P53 26	Unk P53 26	Unk P53 28	Unk P53 28	Unk P53 28	Unk P53 30	Unk P53 30	Unk P53 30	NTC	NTC	NTC
D	Unk EGFR 26	Unk EGFR 26	Unk EGFR 26	Unk EGFR 28	Unk EGFR 28	Unk EGFR 28	Unk EGFR 30	Unk EGFR 30	Unk EGFR 30	NTC	NTC	NTC
E	Unk Actb 29	Unk Actb 29	Unk Actb 29	Unk Actb 29	Unk Actb 29	Unk Actb 29	Unk Actb 30	Unk Actb 30	Unk Actb 30	NTC	NTC	NTC
F	Unk GAPDH 29	Unk GAPDH 29	Unk GAPDH 29	Unk GAPDH 29	Unk GAPDH 29	Unk GAPDH 29	Unk GAPDH 30	Unk GAPDH 30	Unk GAPDH 30	NTC	NTC	NTC
G	Unk P53 27	Unk P53 27	Unk P53 27	Unk P53 29	Unk P53 29	Unk P53 29	Unk P53 30	Unk P53 30	Unk P53 30	NTC	NTC	NTC
H	Unk EGFR 27	Unk EGFR 27	Unk EGFR 27	Unk EGFR 29	Unk EGFR 29	Unk EGFR 29	Unk EGFR 30	Unk EGFR 30	Unk EGFR 30	NTC	NTC	NTC

	1	2	3	4	5	6	7	8	9	10	11	12
A	Unk Actb 14	Unk Actb 14	Unk Actb 14	Unk Actb 16	Unk Actb 16	Unk Actb 16	Unk Actb 18	Unk Actb 18	Unk Actb 18	Unk Actb 18	Unk Actb 18	Unk Actb 18
B	Unk GAPDH 14	Unk GAPDH 14	Unk GAPDH 14	Unk GAPDH 16	Unk GAPDH 16	Unk GAPDH 16	Unk GAPDH 18	Unk GAPDH 18	Unk GAPDH 18	Unk GAPDH 18	Unk GAPDH 18	Unk GAPDH 18
C	Unk P53 14	Unk P53 14	Unk P53 14	Unk P53 16	Unk P53 16	Unk P53 16	Unk P53 18	Unk P53 18	Unk P53 18	Unk P53 18	Unk P53 18	Unk P53 18
D	Unk EGFR 14	Unk EGFR 14	Unk EGFR 14	Unk EGFR 16	Unk EGFR 16	Unk EGFR 16	Unk EGFR 18	Unk EGFR 18	Unk EGFR 18	Unk EGFR 18	Unk EGFR 18	Unk EGFR 18
E	Unk Actb 15	Unk Actb 15	Unk Actb 15	Unk Actb 17	Unk Actb 17	Unk Actb 17	Unk Actb 19	Unk Actb 19	Unk Actb 19	NTC	NTC	NTC
F	Unk GAPDH 15	Unk GAPDH 15	Unk GAPDH 15	Unk GAPDH 17	Unk GAPDH 17	Unk GAPDH 17	Unk GAPDH 19	Unk GAPDH 19	Unk GAPDH 19	NTC	NTC	NTC
G	Unk P53 15	Unk P53 15	Unk P53 15	Unk P53 17	Unk P53 17	Unk P53 17	Unk P53 19	Unk P53 19	Unk P53 19	NTC	NTC	NTC
H	Unk EGFR 15	Unk EGFR 15	Unk EGFR 15	Unk EGFR 17	Unk EGFR 17	Unk EGFR 17	Unk EGFR 19	Unk EGFR 19	Unk EGFR 19	NTC	NTC	NTC

	20	20	20	22	22	22	24	24	24	24	24	24
C	Unk P53 20	Unk P53 20	Unk P53 20	Unk P53 22	Unk P53 22	Unk P53 22	Unk P53 24	Unk P53 24	Unk P53 24	Unk P53 24	Unk P53 24	Unk P53 24
D	Unk EGFR 20	Unk EGFR 20	Unk EGFR 20	Unk EGFR 22	Unk EGFR 22	Unk EGFR 22	Unk EGFR 24	Unk EGFR 24	Unk EGFR 24	Unk EGFR 24	Unk EGFR 24	Unk EGFR 24
E	Unk Actb 21	Unk Actb 21	Unk Actb 21	Unk Actb 23	Unk Actb 23	Unk Actb 23	Unk Actb 25	Unk Actb 25	Unk Actb 25	Unk Actb 25	Unk Actb 25	Unk Actb 25
F	Unk GAPDH 21	Unk GAPDH 21	Unk GAPDH 21	Unk GAPDH 23	Unk GAPDH 23	Unk GAPDH 23	Unk GAPDH 25	Unk GAPDH 25	Unk GAPDH 25	Unk GAPDH 25	Unk GAPDH 25	Unk GAPDH 25
G	Unk P53 21	Unk P53 21	Unk P53 21	Unk P53 23	Unk P53 23	Unk P53 23	Unk P53 25	Unk P53 25	Unk P53 25	Unk P53 25	Unk P53 25	Unk P53 25
H	Unk EGFR 21	Unk EGFR 21	Unk EGFR 21	Unk EGFR 23	Unk EGFR 23	Unk EGFR 23	Unk EGFR 25	Unk EGFR 25	Unk EGFR 25	Unk EGFR 25	Unk EGFR 25	Unk EGFR 25

Good DATA



A new standard for qPCR

Clinical Chemistry 55:4
611–622 (2009)

Special Report

The MIQE Guidelines: *Minimum Information for Publication of Quantitative Real-Time PCR Experiments*

Stephen A. Bustin,^{1*} Vladimir Benes,² Jeremy A. Garson,^{3,4} Jan Hellemans,⁵ Jim Huggett,⁶
Mikael Kubista,^{7,8} Reinhold Mueller,⁹ Tania Nolan,¹⁰ Michael W. Pfaffl,¹¹ Gregory L. Shipley,¹²
Jo Vandesompele,⁵ and Carl T. Wittwer^{13,14}

The MIQE guidelines... “target the **reliability** of results to help ensure the **integrity** of the scientific literature, promote **consistency** between laboratories, and increase experimental **transparency**.”

>4000 Citations

What is MIQE: It's a Checklist

Table 1. MIQE checklist for authors, reviewers, and editors.^a

Item to check	Importance	Item to check	Importance
Experimental design		qPCR oligonucleotides	
Definition of experimental and control groups	E	Primer sequences	E
Number within each group	E	RTPimerDB identification number	D
Assay carried out by the core or investigator's laboratory?	D	Probe sequences	D ^d
Acknowledgment of authors' contributions	D	Location and identity of any modifications	E
Sample		Manufacturer of oligonucleotides	D
Description	E	Purification method	D
Volume/mass of sample processed	D	qPCR protocol	
Microdissection or macrodissection	E	Complete reaction conditions	E
Processing procedure	E	Reaction volume and amount of cDNA/DNA	E
If frozen, how and how quickly?	E	Primer, (probe), Mg ²⁺ , and dNTP concentrations	E
If fixed, with what and how quickly?	E	Polymerase identity and concentration	E
Sample storage conditions and duration (especially for FFPE ^b samples)	E	Buffer/kit identity and manufacturer	E
Nucleic acid extraction		Exact chemical composition of the buffer	D
Procedure and/or instrumentation	E	Additives (SYBR Green I, DMSO, and so forth)	E
Name of kit and details of any modifications	E	Manufacturer of plates/tubes and catalog number	D
Source of additional reagents used	D	Complete thermocycling parameters	E
Details of DNase or RNase treatment	E	Reaction setup (manual/robotic)	D
Contamination assessment (DNA or RNA)	E	Manufacturer of qPCR instrument	E
Nucleic acid quantification	E	qPCR validation	
Instrument and method	E	Evidence of optimization (from gradients)	D
Purity (A ₂₆₀ /A ₂₈₀)	D	Specificity (gel, sequence, melt, or digest)	E
Yield	D	For SYBR Green I, C _q of the NTC	E
RNA integrity: method/instrument	E	Calibration curves with slope and y intercept	E
RIN/RQI or C _q of 3' and 5' transcripts	E	PCR efficiency calculated from slope	E
Electrophoresis traces	D	CIs for PCR efficiency or SE	D
Inhibition testing (C _q dilutions, spike, or other)	E	r ² of calibration curve	E
Reverse transcription		Linear dynamic range	E
Complete reaction conditions	E	C _q variation at LOD	E
Amount of RNA and reaction volume	E	CIs throughout range	D
Priming oligonucleotide (if using GSP) and concentration	E	Evidence for LOD	E
Reverse transcriptase and concentration	E	If multiplex, efficiency and LOD of each assay	E
Temperature and time	E	Data analysis	
Manufacturer of reagents and catalogue numbers	D	qPCR analysis program (source, version)	E
C _q s with and without reverse transcription	D ^c	Method of C _q determination	E
Storage conditions of cDNA	D	Outlier identification and disposition	E
qPCR target information		Results for NTCs	E
Gene symbol	E	Justification of number and choice of reference genes	E
Sequence accession number	E	Description of normalization method	E
Location of amplicon	D	Number and concordance of biological replicates	D
Amplicon length	E	Number and stage (reverse transcription or qPCR) of technical replicates	E
In silico specificity screen (BLAST, and so on)	E	Repeatability (intraassay variation)	E
Pseudogenes, retrotransposons, or other homologs?	D	Reproducibility (interassay variation, CV)	D
Sequence alignment	D	Power analysis	D
Secondary structure analysis of amplicon	D	Statistical methods for results significance	E
Location of each primer by exon or intron (if applicable)	E	Software (source, version)	E
What splice variants are targeted?	E	C _q or raw data submission with RDML	D

^a All essential information (E) must be submitted with the manuscript. Desirable information (D) should be submitted if available. If primers are from RTPimerDB, information on qPCR target, oligonucleotides, protocols, and validation is available from that source.

^b FFPE, formalin-fixed, paraffin-embedded; RIN, RNA integrity number; RQI, RNA quality indicator; GSP, gene-specific priming; dNTP, deoxynucleoside triphosphate.

^c Assessing the absence of DNA with a no-reverse transcription assay is essential when first extracting RNA. Once the sample has been validated as DNA free, inclusion of a no-reverse transcription control is desirable but no longer essential.



^d Disclosure of the probe sequence is highly desirable and strongly encouraged; however, because not all vendors of commercial predesigned assays provide this information, it cannot be an essential requirement. Use of such assays is discouraged.

- qPCR community driven guidelines for essential and desired information in literature;

- Experimental Design
- Sample Information
- Nucleic Acid Extraction
- Reverse Transcription
- qPCR Target Information
- qPCR Oligonucleotides
- qPCR Protocol
- qPCR Validation
- Data Analysis

Temperature Gradient

dynamic thermal gradient



Test 12 sets of primer
Life Science Group

- MIQE guideline
–Evidence of optimization

Table 1. MIQE checklist for authors, reviewers, and editors.^a

Item to check	Importance	Item to check	Importance
Experimental design		qPCR oligonucleotides	
Definition of experimental and control groups	E	Primer sequences	E
Number within each group	E	RTPimerID: Identification number	D
Assay carried out by the core or investigator's laboratory?	D	Probe sequences	D ²
Acknowledgment of authors' contributions	D	Location and identity of any modifications	E
Sample		Manufacturer of oligonucleotides	
Description	E	Purification method	D
Volume/mass of sample processed	D	qPCR protocol	
Microdissection or macrodissection	E	Complete reaction conditions	E
Processing procedure	E	Reaction volume and amount of cDNA/DNA	E
if frozen, how and how quickly?	E	Primer, (probe), Mg ²⁺ , and dNTP concentrations	E
if fixed, with what and how quickly?	E	Polymerase identity and concentration	E
Sample storage conditions and duration (especially for FFPE ³ samples)	E	Buffer/kits identity and manufacturer	E
Nucleic acid extraction		Exact chemical composition of the buffer	
Procedure and/or instrumentation	E	Additives (SYBR Green I, DMSO, and so forth)	E
Name of kit and details of any modifications	E	Manufacturer of plates/tubes and catalog number	D
Source of additional reagents used	D	Complete thermocycling parameters	E
Details of DNase or RNase treatment	E	Reaction setup (manual/robotic)	D
Contamination assessment (DNA or RNA)	E	Manufacturer of qPCR instrument	E
Nucleic acid quantification		qPCR validation	
Instrument and method	E	Evidence of optimization (from gradients)	D
Purity (A ₂₆₀ /A ₂₃₀)	D	Specificity (gel, sequence, melt, or digest)	E
Yield	D	For SYBR Green I, C _q of the NTC	E
RNA integrity: method/instrument	E	Calibration curves with slope and y intercept	E
RIN/ROU or C _q of 3' and 5' transcripts	E	PCR efficiency calculated from slope	E
Electrochromic traces	D	C _q for PCR efficiency or SE	D

qPCR validation	
Evidence of optimization (from gradients)	D
Specificity (gel, sequence, melt, or digest)	E

Temperature and time	E	Data analysis	
Manufacturer of reagents and catalogue numbers	D ²	qPCR analysis program (source, version)	E
C _q with and without reverse transcription	D ²	Method of C _q determination	E
Storage conditions of cDNA	D	Outlier identification and disposition	E
qPCR target information		Results for NTCs	
Gene symbol	E	Justification of number and choice of reference genes	E
Sequence accession number	E	Description of normalization method	E
Location of amplicon	D	Number and concordance of biological replicates	D
Amplicon length	E	Number and stage (reverse transcription or qPCR) of technical replicates	E
In silico specificity screen (BLAST, and so on)	E	Repeatability (intraassay variation)	E
Pseudogenes, retrospseudogenes, or other homologs?	D	Reproducibility (interassay variation, CV)	D
Sequence alignment	D	Power analysis	D
Secondary structure analysis of amplicon	D	Statistical methods for results significance	E
Location of each primer by exon or intron (if applicable)	E	Software (source, version)	E
What splice variants are targeted?	E	C _q or raw data submission with RDML	D

Bio-Rad's CFX series Real-time PCR Family

Price



CFX Connect

- Mid-price
- 96-well
- Fast block
- 2 targets



CFX96 Touch

- Premium
- 96-well
- Fast block
- 5 targets
- Stand-alone run




CFX384 Touch


- Premium
- 384-well
- 4 targets
- Stand-alone
- High throughput
- Stand-alone run

Features & Flexibility


CFX Connect™



- 3 filtered LEDs
- 3 filtered photodiodes
- 2 targets multiplex




- Temperature Gradient
- High resolution melt





- 96 well 1-50 μ L
- Max ramp rate 5°C/s






- Optics shuttle transverses the plate at the completion of each data acquisition step
- Optics shuttle centered above each well resulting in identical light path and intensity for each well
- No need to sacrifice one channel to data collection of a passive reference dye





Filtered LEDs



Filtered photodiodes



Strength of our optic module

Accurate data collection


Consistent light source

Minimal maintenance

Data collected from every well


Patented Block Design

- Accuracy: $\pm 0.2^{\circ}\text{C}$
- Uniformity: $\pm 0.4^{\circ}\text{C}$ within 10 seconds



Accuracy

Programmed temperature
= Actual well temperature





Uniformity

All wells in block are
same temperature


Performance Thermal Cycler Gradient

- Annealing temp is critical for reaction specificity and efficiency
- Save time and reagents optimizing annealing temp in one run
 - Program up to a 24°C gradient, back-to-front (back row is hotter)
 - 30-100°C range
 - “Dynamic Ramping” – Wells hit set temperature point together, and then maintain the same dwell to eliminate time as a reaction variable



columns on C1000 with 96Fast reaction block


Assay optimization




	For 1	Rev 1	For 2	Rev 2	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
A	70.0	U	U					U	U																				
B	69.7	U	U					U	U																				
C	69.2	U	U					U	U																				
D	68.6	U	U					U	U																				
E	67.5	U	U					U	U																				
F	66.5	U	U					U	U																				
G	65.0	U	U					U	U																				
H	63.5	U	U					U	U																				
I	61.7	U	U					U	U																				
J	60.1	U	U					U	U																				
K	58.8	U	U					U	U																				
L	57.6	U	U					U	U																				
M	56.6	U	U					U	U																				
N	55.9	U	U					U	U																				
O	55.3	U	U					U	U																				
P	55.0	U	U					U	U																				

10° above design (rows A-F)



5° below design (rows G-P)



- 70.0
- 69.7
- 69.2
- 68.6
- 67.5
- 66.5
- 65.0
- 63.5
- 61.7
- 60.1
- 58.8
- 57.6
- 56.6
- 55.9
- 55.3
- 55.0




Optimal Annealing Range




Overcome non-specific amplification

	1	2	3	4
A	68.0		Unk	Unk
B	67.0		Unk	Unk
C	64.9		Unk	Unk
D	61.4		Unk	Unk
E	57.0		Unk	Unk
F	53.5		Unk	Unk
G	51.2		Unk	Unk
H	50.0		Unk	Unk




Overcome non-specific amplification

	1	2	3	4
A	68.0		Unk	Unk
B	67.0		Unk	Unk
C	64.9		Unk	Unk
D	61.4		Unk	Unk
E	57.0		Unk	Unk
F	53.5		Unk	Unk
G	51.2		Unk	Unk
H	50.0		Unk	Unk



Overcome non-specific amplification



	1	2	3	4
A	68.0		Unk	Unk
B	67.0		Unk	Unk
C	64.9		Unk	Unk
D	61.4		Unk	Unk
E	57.0		Unk	Unk
F	53.5		Unk	Unk
G	51.2		Unk	Unk
H	50.0		Unk	Unk





Gradient Optimization of Annealing Temperature

- Serial dilutions run at 8 temps, ranging from 55°C to 68°C
- Reactions at 62°C annealing have low Cts and highest reaction efficiency



67°C




62°C



56°C





Use of Synthetic Templates

Synthetic templates allow for

- Sample independent validation of assay performance
- Comparison across plates (inter-run calibrators)
- Compare across labs

*Nucleic Acids Research, 2009, 1-9
doi:10.1093/nar/gkp721*


External oligonucleotide standards enable cross laboratory comparison and exchange of real-time quantitative PCR data

Joëlle Vermeulen¹, Filip Pattyn¹, Katleen De Preter¹, Liesbeth Vercruysse¹,
Stefaan Derveaux¹, Pieter Mestdagh¹, Steve Lefever¹, Jan Hellemans^{1,2},
Frank Speleman¹ and Jo Vandesompele^{1,2,*}

RNA quality


- RNA quality analysis methodology
 - Analysis A260/A280 ratio and rRNA bands on agarose gel.
 - Automated microfluidic capillary electrophoresis system.
 - 3`5` assay
 - SPUD assay (absence of PCR inhibitors)
 - NRT PCR assay (absence of DNA contamination)

- Analyze RNA quantity and quality
- UV absorbance
 - $A_{260}/A_{280}=2.0$,
 - A_{260} : 0.15~1
 - $A_{260}/A_{230}=2.5$
 - Gel electrophoresis analysis
28S rRNA, 18S rRNA



28S:18S≈2:1

RNA quality




SPUD: detection presence of inhibitors

RNA Quality Control Assay

RNA Quality Assay – Two assays (RQ1 and RQ2) that target the same transcript, but at two different locations. A significantly higher presence of the shorter transcript (RQ1) indicates the RNA is degraded. A difference in C_q less than 3 between the two assays indicates acceptable RNA quality.

RQ1




PASS

$$= \Delta C_q \leq 3$$

RNA degradation will not negatively affect results!

RQ2



FAIL

$$= \Delta C_q > 3$$

RNA is degraded and could be negatively affecting results – most likely need to redo sample extraction

Bio-Rad RT Kit





All Bio-Rad's iScript Kit

87

All Bio-Rad's iScript Kits:

- Use an MMLV enzyme with RNase H⁺ activity
- Come with an optimal blend of oligo dTs and random primers, (except the *Select* kit, which provides the primers in a separate tube)
- Include a potent RNase inhibitor to protect the RNA from degradation prior to cDNA synthesis
- Were benchmarked against the top selling competitor products to demonstrate superior sensitivity and dynamic range








Review – Best Product for the Job


	iScript Reverse Transcription Supermix for RT-qPCR	iScript gDNA Clear cDNA Synthesis Kit	iScript Advanced cDNA Synthesis Kit for RT-qPCR	iScript Select cDNA Synthesis Kit	iScript cDNA Synthesis Kit
Features	Minimal reaction setup time	Effective gDNA removal before reverse transcription	Maximum RNA input for high cDNA yields	Flexible priming options	Reliable value solution
Cost	\$\$	\$\$\$	\$\$\$	\$	\$
Total input RNA Range	1 µg - 1pg	1 µg - 1pg	7.5 µg – 100 fg	1 µg - 1pg	1 µg - 100 fg
Kit Format	1 tube	3 tubes	2 tubes	5 tubes	2 tubes
Kit Contents	<ul style="list-style-type: none"> • 5x iScript RT Supermix • No-RT control supermix 	<ul style="list-style-type: none"> • 5x iScript RT Supermix • No-RT control supermix • DNase • DNase buffer 	<ul style="list-style-type: none"> • iScript Reverse Transcriptase • 5x iScript Advanced Reaction Mix 	<ul style="list-style-type: none"> • iScript Reverse Transcriptase • 5x iScript Reaction Mix • 3 priming options 	<ul style="list-style-type: none"> • iScript Reverse Transcriptase • 5x iScript Reaction Mix
Time to Produce cDNA	26 min	36 min	21 min	40-90 min	26 min

qPCR Supermix




SsoAdvanced™ Universal Supermixes			iTaq™ Universal Supermixes	
				
SsoAdvanced™ Universal SYBR® Green Supermix	SsoAdvanced Universal Probes Supermix	SsoAdvanced Universal Inhibitor-Tolerant Supermix	iTaq™ Universal SYBR® Green Supermix	iTaq Universal Probes Supermix
Robust performance with demanding templates and under challenging conditions			Advanced formulation for reliable, reproducible gene expression analysis on any qPCR instrument	





型號	敘述
Real-Time PCR 耗材	
八連排	
TCS-0803	Optical Flat 8-Cap Strips, 960 reactions
TLS-0851	Low-Profile 0.2 ml 8-Tube Strips without Caps, 960 reactions
TLS-0801	Low-Profile 0.2 ml 8-Tube Strips without Caps, 960 reactions
96 或 48 孔盤	
MSB-1001	Microseal 'B' Adhesive Seals, Pkg of 100
MLL-9651	Multiplate White Low-Profile 96-Well Unskirted PCR Plates, Pkg of 25
MLL-9601	Multiplate Low-Profile 96-Well Unskirted PCR Plates, Pkg of 25
HSP-9955	Hard-Shell Thin-Wall 96-Well Skirted PCR Plates, Pkg of 50
MLL-4851	Multiplate White Low-Profile 48-Well Unskirted PCR Plates, Pkg of 50



反轉錄相關試劑

170-8890	iScript™ cDNA Synthesis Kit, 25 reactions
170-8891	iScript™ cDNA Synthesis Kit, 100 reactions
170-8896	iScript™ Select cDNA Synthesis Kit, 25 reactions
170-8897	iScript™ Select cDNA Synthesis Kit, 100 reactions
170-8892	iScript™ One-Step RT-PCR Kit With SYBR® Green, 50 reactions
170-8893	iScript™ One-Step RT-PCR Kit With SYBR® Green, 200 reactions
170-8894	iScript™ One-Step RT-PCR Kit for Probes, 50 reactions
170-8895	iScript™ One-Step RT-PCR Kit for Probes, 200 reactions

Real-Time PCR 混合液

172-5848	iQ™ Multiplex Powermix, 50 reactions
172-5849	iQ™ Multiplex Powermix, 200 reactions
170-8860	iQ™ Supermix, 100 reactions
170-8862	iQ™ Supermix, 500 reactions
170-8864	iQ™ Supermix, 1,000 reactions
170-8880	iQ™ SYBR® Green Supermix, 100 reactions
170-8882	iQ™ SYBR® Green Supermix, 500 reactions
170-8884	iQ™ SYBR® Green Supermix, 1,000 reactions
170-8885	iQ™ SYBR® Green Supermix, 2,000 reactions
172-5200	SsoFast™ EvaGreen® Supermix, 200 reactions
172-5201	SsoFast™ EvaGreen® Supermix, 500 reactions
172-5202	SsoFast™ EvaGreen® Supermix, 1,000 reactions
172-5203	SsoFast™ EvaGreen® Supermix, 2,000 reactions
172-5230	SsoFast Probes Supermix, 200 reactions
172-5231	SsoFast Probes Supermix, 500 reactions
172-5232	SsoFast Probes Supermix, 1,000 reactions


Industry-leading CFX Manager Software




▪ User-friendly




▪ Effective Wizard




▪ Concise Interface



Data Collecting and Analysis



Powerful Function modules



▪ Rigorous algorithm

CFX Manager software





- Absolute quantification
- Relative quantification (Gene expression analysis)
- Melt curve analysis
- Genotyping
- QC
- Run information
- Qualification Plate Run




CFX Manager software

- Sample group
- PCR efficiency correction
- Normalization with multiple reference genes
- Gene Study (combination of multiple plates)
- Screening best reference genes by M value
- Inter-run-calibration
- Comply with MIQE guideline
 - RDML export RDML
- Visualize your data
 - Bar chart
 - Cluster gram
 - Scatter plot
 - Volcano plot
 - Heat Map
 - Results table





Test primers of both target gene(s) and reference gene(s)





For each gene involved respectively

Inter-run calibrators


Correct analysis model





Challenge of q-PCR





Amplification



What about an new method - Digital PCR




PCR reaction that is partitioned.




Many thousands
of discrete measurements


ddPCR Principle



- End point (0's or 1's)
- Less sensitive to PCR efficiency
- No standard curve
- More tolerant to PCR inhibitors



The sample is partitioned into 20,000 droplets, with target and background DNA randomly distributed among the droplets.




After PCR amplification, each droplet provides a fluorescent positive or negative signal indicating the target DNA was present or not present after partitioning. Each droplet provides an independent digital measurement.

“X” target copies

Positive and negative droplets are counted for the sample and the software calculates the concentration of target DNA as copies per microliter.

Workflow of Droplet Digital PCR


1. Make Droplets



Droplet Generator




2. Cycle Droplets



Bulk PCR
Thermal Cycler



3. Read Droplets



Droplet Reader





“X” target
copies

Results


- Readout: copies/ μ l
- Dynamic range: 1–100,000 copies/well
(~330 ng human genomic DNA)

ddPCR assays are more easy to design and validate




- More Tolerance to:
1. Non-specific binding
 2. Primer dimer
 3. PCR inhibitor

Non-specific cross reaction




qPCR



ddPCR

Temperature Gradient: EvaGreen


- Primer dimers can be visualized with EvaGreen chemistry
- Primer dimer frequency increases at lower annealing temperature and are visible in NTCs




- More Tolerance to:
1. Non-specific binding
 2. Primer dimer
 3. PCR inhibitor

Applications

Copy Number Variation




Mutation Detection



Stem Cell Analysis


Virus Detection



NGS Cross Validation



Gene Expression
mRNA, small RNA, lncRNA



Q&A



IT STARTS WITH A
DROPLET
AND ENDS IN
DISCOVERY.
Droplet Digital™ PCR

E-mail: charlene_lu@bio-rad.com
Phone: 0987-792-623