


# 結核菌（痰液）培養採集須知

1.採檢前	2.採檢過程	3.送檢流程
 <p>圖片取自 PIXTA 網站</p>  <p><a href="https://www.google.com.tw/search?safe=off&amp;rlz=1C1AOHY_zh-TWTW708TW708&amp;biw=1280&amp;bih=889&amp;tbn=isch&amp;sa=1&amp;q=%E7%97%B0&amp;oeq=%E7%97%B0&amp;gs_l=img.3..0110.2413.3880.0.5053.2.2.0.0.0.66.112.2.0..0..1..1.4.64.img..0.2.110.4DVXW1N1WV..0.0.0.0">https://www.google.com.tw/search?safe=off&amp;rlz=1C1AOHY_zh-TWTW708TW708&amp;biw=1280&amp;bih=889&amp;tbn=isch&amp;sa=1&amp;q=%E7%97%B0&amp;oeq=%E7%97%B0&amp;gs_l=img.3..0110.2413.3880.0.5053.2.2.0.0.0.66.112.2.0..0..1..1.4.64.img..0.2.110.4DVXW1N1WV..0.0.0.0</a></p>	 <p><a href="https://www.google.com.tw/search?safe=off&amp;rlz=1C1AOHY_zh-TWTW708TW708&amp;biw=1280&amp;bih=889&amp;tbn=isch&amp;sa=1&amp;q=sputum&amp;oeq">https://www.google.com.tw/search?safe=off&amp;rlz=1C1AOHY_zh-TWTW708TW708&amp;biw=1280&amp;bih=889&amp;tbn=isch&amp;sa=1&amp;q=sputum&amp;oeq</a></p> 	 <p style="text-align: center;">↓ 盡 快 送 驗 ↓</p> <p style="text-align: center;">病理檢驗科</p> 
<p><b>注意事項：</b></p> <p>①以清晨第一口痰為佳。</p> <p>②檢體不可有菜渣、口水污染。</p>	<p>①先以開水漱口 ↓</p> <p>②用力咳出肺部深處 ↓</p> <p>③咳入無菌容器內 ↓</p> <p>④將螺旋蓋蓋好 ↓</p> <p>⑤放入塑膠夾鏈袋</p>	<p>⑥未能一小時內送驗，需放置在 2-8℃ 冰箱冷藏。 *結核菌檢查可冷藏 &lt; 72 小時，但請採檢第三天當天送檢。</p> <p style="text-align: center;">↓</p> <p>⑦請當日與檢查單送至病理檢驗科。</p>

如有疑問請洽病理檢驗科諮詢，諮詢電話:08-7557885 分機84830、84840

屏東榮民總醫院永遠關心您